


LVL'D UP

The next Generation



Redesigned 2D tubes

- 1** 2D bottom code / Protective rim against mechanical influences (scratches)
- 2** 1D side code / HR side code
- 3** Open space for customized coding
- 4** SAFE® low leachables medical polypropylene



New cap design

Excellent sealing
Easy manual handling
Automation-friendly cap



Your free sample pack


Use the QR code or email us at
info@lvl-technologies.com



Experience our new LVL Racks




3 new areas for additional rack and tube information



Easy Lock

One simple move to open and close



Anti-freeze gradient

Prevents ice formation



Tube Centering

Automation-friendly beveled domes



Curved well rim

Safety guides the tube and reduces collisions