



CAP THEM ALL!

Now available! LVL semi-automated hand-held decapper / capper. The best choice for medium throughput sample storage workflows.

- Available in SBS 96, 48 and 24 formats for all models of LVL SAFE® 2D barcoded tubes
- Rapid capping and decapping
- Ejects caps safely



Spring loaded retention balls keep caps securely in place.



Technical specifications located on back

Find out more:
LVL cappers online




SAFE®
Store Aliquots For Ever

Order now.

lvl-technologies.com/en/devices/capper--decapper

Optional accessory LVL Cap Rack



Conveniently store tube caps on our cap racks.

SBS 96: RCK-CR96-P-L
SBS 48: RCK-CR48-P-L
SBS 24: RCK-CR24-P-L



4-channel
semi-automatic capper



6-channel
semi-automatic capper



8-channel
semi-automatic capper

Rack format	SAFE® SBS 24	SAFE® SBS 48	SAFE® SBS 96
Opening / closing speed	0.5s / 4.7s	0.5s / 4s	0.5s / 3.5s
Motors	4 individually controlled high precision motors *	6 individually controlled high precision motors *	8 individually controlled high precision motors *
Stand included	✓	✓	✓
Article number	CDC-4CH-XT	CDC-6CH-XT	External thread: CDC-8CH-XT-4 Internal thread: CDC-8CH-IT-2

* 4 - 8 individual motors (depending on the model) and standardized closing torques enable reliable and reproducible closing results.