


# The LVL Sample Tubes

For the daily use of sample handling and archiving

- Material: translucent polypropylene (PP)
- Available as individual Single Tubes or 8-Tube Strips for use with 96 SBS Tube Racks
- Volume 0.65 ml or 1.2 ml with round bottom
- Very low protein binding therefore well-suited for sample dilution
- Closable either with a cover mat, 8-Tube Strip or self-adhesive foil
- Autoclavable for 20 min at 121 °C
- Temperatur-restistant to -80 °C



## Tubes

Article No.	Volume	Quantity per packaging unit
226.01.BL	1,2 ml	5.000 pcs
223.01.BL	0,65 ml	5.000 pcs
226.08.BL	1,2 ml	600 pcs
223.08.BL	0,65 ml	600 pcs

1 Single Tubes 0,65 ml and 1,2 ml • 2 8-Tube Strips 1,2 ml and 0,65 ml • 3 Normal Rack • 4 High Rack