

SAFE® ACCESS Tube Picker

Semi-automated tube picker



FEATURES

- Fast and reliable tube picking from 96 position racks
- Works with frozen and thawed tubes
- Create picklists in excel, csv, text, xml or json
- A 2D barcode scanner for tube confirmation and for generating picklists directly from the tube barcodes
- 1D scanner for rack barcode included as standard
- DataPaq™ Software with free lifetime upgrades

Easy to Set Up and Easy to Use

The SAFE® Access Tube Picker works out of the box and needs no set up or calibration. Instrument control is through user friendly software which provides an intuitive interface for easy operation.

Tubes can be picked from a single rack or a picklist can be set up to select tubes across multiple racks. Pick list operations can process racks from the linear barcode or manually select wells to pick.

SPECIFICATIONS

- Dimensions = 300mm x 250mm x 180mm (WxDxH)
- Total Weight = 5.7kg
- Power adaptor input = 100-240VAC
- Power adaptor output = 12VDC, 5A
- Picklist import file formats are excel, text, csv, json and xml
- Operating system on host PC = Windows 7, 8, 10
- Compatible with 96 position SBS rack

INCLUDED WITH THE READER

- SAFE® Access Tube Picker instrument
- Power adaptor, power cable
- Software and License
- 1D rack code reader and bracket
- User guide

OPERATION

AUTOMATIC MODE

Designed to allow for operation of the Tube Picker with minimal interaction on control PC.

Worklist File

A file that describes the tubes that are to be selected, there are multiple formats available.

Ordering Information

Article No.	Description	Quantity
TPI96	SAFE® ACCESS Tube Picker	1