

# SAFE® READ Single Tube

## SMALL AND COST EFFECTIVE 2D BARCODE READER



## FEATURES

- Fast scanning of 2D tubes
- Compatible with all 2D Datamatrix tubes
- 'Keyboard Wedge' setting allows data entry directly into your application
- Strong mineral glass window
- Software with free lifetime upgrades
- Cryoprotection option

## Easy to Set Up and Easy to Use

The SAFE® READ Single Tube reader offers easy, out-of-the-box set up and is delivered pre-calibrated and ready to read.

This instrument can read all commercially available 2D datamatrix tubes.

The SAFE® READ Single Tube reader can be used for applications that require data input without direct interaction with your computer.

# SPECIFICATIONS

- Dimensions = 60mm x 90mm x 60mm (WxDxH)
- Total Weight = 0.4kg
- Powered by USB - connect directly into computer
- Operating system on host PC = Windows 7, 8, 10
- Warranty period = 2 years
- Free Software upgrades for life
- Software includes remote control capability for integration applications

# INCLUDED WITH THE READER

- SAFE® READ Single Tube instrument
- Power adaptor and country relevant power cable
- Software
- User guide

# OPTIONS

## CRYOPROTECTION (CP)

Prevents the reader window from misting up when a cold rack from a freezer is scanned.  
Our passive technology does not use heat or blown air.

## 1D LINEAR SCANNER

Reads the 1D code on the edge of the rack holding the tubes

# Ordering Information

Article No.	Description	Quantity
DMST	SAFE® READ Single Tube	1
DMST-CP	SAFE® READ Single Tube with Cryoprotection	1
BCRR	SAFE® 1D Rack scanner (in combo with 2D Rack Reader)	1